

Device Specifications	VIAFLO 96	VIAFLO 384
Dimensions (W x D x H)	42 cm x 30 cm x 54 cm (17" x 12" x 21")	
Weight	25.7 kg (56.7 lbs)	27.7 kg (61.1 lbs)
Power requirements	100 - 240 VAC, 50/60 Hz	
Operating conditions	5 °C - 40 °C, <80 % Humidity, non-condensing	
Pipetting channels	24, 96	24, 96, 384
Pipetting technology	air displacement	
Compatible formats	24-, 96-, and 384-well microplates, shallow and deep well	24-, 96-, 384, and 1536-well microplates, shallow and deep well