

Device Specifications	VIAFLO 96	VIAFLO 384
Dimensions (W x D x H)	42 cm x 30 cm x 54 cm (17" x 12" x 21")	
Weight	25.7 kg (56.7 lbs)	27.7 kg (61.1 lbs)
Power requirements	100 – 240 VAC, 50/60 Hz	
Operating conditions	5 °C – 40 °C, <80 % Humidity, non-condensing	
Pipetting channels	96	384
Pipetting technology	air displacement	
Compatible formats	96- and 384-well microplates, shallow and deep well	